

Control Unit CP5000

Compatible with
ATX3000 product range.

With its 15" TFT flat screen the CP5000 can display text addresses for call types, locations, resident/patient names, active call list and alarm location with an optional floor plan/map.

The status of calls is displayed and more information can be requested where necessary, for example, brief patient medical notes can be called up.

Five alarm tones are used; call, assist, emergency, reminder and maintenance to make the system as user friendly and simple to operate. Each of these alarm tones can be altered on programming thus future proofing the system to any changes in alarm tone regulation.

'Reminder' call's can be staff programmed to automatically activate at preset intervals for use in the monitoring of nightly checks or to assist with timed medication.

Multiple zones can be programmed and given names on each panel. The unit has an integral call logger giving accurate information of all patient or staff activated events as well as maintenance information. This information can also be outputted to our 'Call Logger' allowing greater manipulation of the information.

There is a power-fail indicator and an internal battery provides around 4 hours back up power.

Interfaces are provided for paging, additional Base Receivers, Mimic panels and Audio Visual Indicators (over door lights).



Aid Call
HEALTHCARE COMMUNICATIONS

Control Unit CP5000 System Features

Aid Call Ltd Healthcare Communications • Linhay House • Linhay Business Park • Ashburton • Devon TQ13 7UP
FREEPHONE 0800 052 3616 • T 01364 651593 • F 01364 651451 • E healthcaresales@aidcall.co.uk
www.aidcall.co.uk

Operational Modes

- Nurse call
- Warden call

Display

- 15 inch TFT display with 120 degree viewing angle
- Text address for call location and resident names
- Patient notes/reminders
- Active call list and staff location list
- Mapping/floor plan (optional)
- Bespoke welcome screen
- Mains power failure indication

Sounder

- Distinct alarm tones for call, assist, emergency, reminder and maintenance
- Optional Alarm tones for all alarm types
- Adjustable volume control
- Call mute facility (programmable)

Function Buttons

- Depending on the operation mode the buttons are:
Volume; Mute; Call accept; ; Map; Resident/patient Info;
- A Cursor button to provide basic navigation of menu driven features

Call History

- Inbuilt Call Logger
- Optional Interface to Aid Call PC call logger application

Paging

Paging information consists of zone message, alarm transmitter indent and call type

- Interface for multiple paging protocols
- Programmable re-page delay 15-600 seconds
- Paging programmable by zone
- Group codes for: Call, Assistance, Emergency, Night group, Maintenance

Base Receivers

Each RX3000 Base Receiver has a unique ident, any failure is identified on the screen

- Up to 4 Base Receivers can be connected

Active Call Upgrade

- Call tone incremented to next level if call remains unanswered for 2 minutes (Optional)

Common Call Type

- Multiple zoning capabilities
- Multiple common call facility
- Multiple night setting facility
- Multiple call type which will be displayed on selected zones,
e.g. emergency, assistance, call, maintenance, reminder and nurse presence

Maintenance Notification

- Indication if a transmitter is not working.
- Low battery indication for individual transmitters

Mimic Displays

- Hard wired and radio mimic panels may be linked to the system to provide call indication anywhere in the building

Dimensions

Approx 40cm by 34cm by 11cm

Weight

Approx 4.1kg without battery
Approx 6.6kg with battery

Technical Information

Frequency 173.225MHz
Approval ETS 300 220

Safety

BS EN 60950 • EMC • R&TTE Directive 1999/5/EC • ETS 301 489-1 • ETS 301-489-3 • RoHS/ WEEE compliant