

Sensatec® Fall Management Solutions



Where falls are a possibility, Sensatec® stands guard, delivering the security and peace of mind everyone deserves.

Sensatec® Fall Management Provides:

“PLUG AND PLAY” SYSTEM SETUP

Enables you to plug in the Sensatec® Sensor Pad, position the patient on the pad, and the system automatically begins to monitor.

ONE PORTABLE UNIT WITH AN EASY TO USE, ERGONOMIC DESIGN

Can be used with a wheelchair or bed.

SMART SYSTEM TECHNOLOGY

Provides audible and visual notification of the current status of the unit.

ENHANCED ALARM NOTIFICATION

The optional nurse call adaptor links to your existing nurse call system for alarm notification.

AUTO RESET

The Sensatec® Fall Management System automatically resets itself after an alarm when the patient's weight is reapplied to the Sensor Pad.

DESCRIPTION	BED SENSOR PADS	CHAIR SENSOR PADS
PART NUMBER	1000-1879K	1000-1869K
7-DAY PADS	(10-PACK)	(10-PACK)
PART NUMBER	1000-1819K	1000-1829K
30-DAY PADS	(10-PACK)	(10-PACK)
PART NUMBER	1000-1849K	1000-1859K
90-DAY PADS	(10-PACK)	(10-PACK)
PART NUMBER	1000-1899K	1000-1889K
180-DAY PADS	(10-PACK)	(10-PACK)
SIZE	30" x 12" x 0.187"	12" x 12" x 0.187"
COVER MATERIAL	POLYETHYLENE	POLYETHYLENE
CORD LENGTH	54"	24"
WARRANTY	7-DAY, 30-DAY, 90-DAY, OR 180-DAY	7-DAY, 30-DAY, 90-DAY, OR 180-DAY

FEATURE/UNIT	ST260	ST270	ST320
PART NUMBER	1000-0875	1000-0876	1000-0878
SILENCE/OFF/ON	•	•	•
PROTECTIVE BOOT	•	•	•
LOW BATTERY NOTIFICATION	•	•	•
AC POWER SUPPLY		•	•
AUTO RESET	•	•	•
AUTO ON	•	•	•
NURSE CALL		•	•
DELAY/VOLUME			•



Visit www.rft.com or contact us directly to purchase Sensatec® products.

RFT TECHNOLOGIES™

3125 N. 126th Brookfield, WI 53005
 800-669-9946 / Ph: 262 790 1771 / Fax: 262.790.1784
www.rft.com / Email: info@rft.com