

Accessories

Additional trigger devices can be added to the Nurse Call system to cater for scenarios such as wandering residents, people with disabilities and particularly those residents who may need additional care and attention.

These accessories are the extra eyes and ears to the system and staff and can include the following additional options:



one Pendant APT3000

one

Pendant APT3000

A radio nurse call transmitter designed to be worn as a pendant.

The pendant incorporates a visual reassurance indicator, large call button and neck cord with break point.

Specification

- Dimensions: 100(L), 45(W), 20(D) mm
- Weight: 85g
- Call Types: 5
- Power: 4 x N size alkaline batteries
- Operating temperature: -12 to +55 degrees Celsius
- Frequency: 173 Mhz
- Specification: MPT1344



two Audio Visual Indicator AVI3000

two

Audio Visual Indicator AVI3000

A multi coloured over door indicator as additional visual indication to calls being activated.

Specification

- Dimensions: 79(L), 86(W), 90(D) mm
- Weight: 85g
- Power: 6-12vDC



Accessories

Aid Call Ltd Healthcare Communications • Linhay House • Linhay Business Park • Ashburton • Devon TQ13 7UP
FREEPHONE 0800 052 3616 • T 01364 651593 • F 01364 651451 • E healthcaresales@aidcall.co.uk
www.aidcall.co.uk

Bed Weight Sensor

In-bed weight sensor monitors change in state when resident leaves bed immediately or after delayed period.

Chair Weight Sensor

In-chair weight sensor monitors change in state when resident leaves chair immediately or after delayed period.



Door Monitoring

Door monitoring allows staff to be alerted when internal or external doors are opened or closed. We are also able to provide the ability to monitor a door bell both with optional override.

Air Switch

Waterproof device to be used to make calls in wet environments i.e. Bath or Shower.

PIR Movement Sensor

Switchable movement detector to protect restricted areas or detect wandering residents.



Ring Detector

To give an indication on the Nurse Call System that a telephone is ringing (more than one can be used per system).

Additional Monitoring Devices

Examples of where additional monitoring devices can be used:

- Epilepsy
- Incontinence
- Environmental controls
- Sip/Puff switches
- Sound activated switches
- Plus a selection of innovative products for people with special needs.